



Read more patterns precisely and quickly. *The next Generation of* **HEp-2**

PRECISE, EFFICIENT, PROVEN


ANA HEp-2 IFA Patterns




ANA HEp2 Homogenous IgG Positive
(400 X)




ANA HEp2 MA IgG Positive
(400 X)




ANA HEp2 MSA IgG Positive
(400 X)




ANA HEp2 Nucleolar IgG Positive
(400 X)




ANA HEp2 Speckled IgG Positive
(400 X)




ANA HEp2 Centromere IgG Positive
(400 X)



ANA HEp2 IgG Negative
(400 X)



ANA HEp2 SSA IgG Positive
(400 X)



ANA HEp2 SSB IgG Positive
(400 X)

Scimedx has created an all new 3rd generation ANA HEp-2 Immunofluorescent test line. These HEp-2 substrate slides offer superior clarified patterns, tremendous sensitivity, and minimal background.

Allowing you to read more patterns precisely and quickly.

- Available in 5 wells plus 2 control wells, or 12 wells plus 2 control wells
- Test sizes 60 (1060L); 120 (1120L); 300 (H300L)
- All components can be sold separately
- Automation ready format

SCIMEDX
CORPORATION



www.scimedx.com